

## WOBURN LOW PROFILING BED



**Product Code:**

HLB797.03

**Description:**

The dependable 4 section profiling bed now features high quality control system, a reduced minimum platform height and a higher load capacity. By lowering the platform to only 220mm from floor level, the risk of injury for those with a tendency to roll out of bed is significantly reduced. Comfort setting of height, back rest and knee break angles is easy when using the lockable handset.

**Safe Working Load:**

220kg (34 Stone)

**Max. Patient Weight:**

185kg (29 Stone)

**Mattress Base Size:**

2000 x 900mm

**Max. Mattress Depth:**

205mm

**Adjustable Height:**

220 - 650mm

**Angle of Back Rest:**

Up to 70 degrees

**Angle of Knee Break:**

Up to 30 degrees

**Bed Weight:**

89.5kg (without side rails)

**Key Information:**

4 Section Profiling Bed  
Electronically Powered Back Rest & Thigh Rest  
Manual Leg Rest Adjustment  
Wraparound Wooden Head and Foot Boards  
Side Rails - Optional  
Trendelenburg Function and Auto-Regression - Optional  
In-Built Extension - Optional  
Lockable Handset  
Low Voltage – 29v; 9v Battery Backup  
Comprehensive Manual  
CE marked

**Supplier Contact Details:**

t: +44(0)1709 377172 e: enquiries@harvesthealthcare.co.uk w: www.harvesthealthcare.co.uk  
Sheaf house, Bradmarsh Way, Bradmarsh Business Park, Rotherham UK, S60

<UJ Ygh<YUTHVWYF @K gi dd'JgYei Jda YbhUxb Udd JubbWg hC Wj YFU Rubj Y cZa YXJW VebXfHcbg' GdYUWwKcJW cZYei Jda Ybhjg hY gc Y YgdcbgMjMmChZkY WbJMU dFUVWfjcbYcf cFUVYF  
cj YGgYyYj' JbXj Jg U dUrhYbVWYF' Rig\_ UggYgga YbhjgdYUWm hC gYFa i ghVY WUfYfj cJ hdrf: hC hYI j gY cZU drcJk Wb'  
PrCxi VbUxb XUa Vebhbi Umi bXy bC cXj Ycda YbhUxb a Umh.YfYZc:fY XjJ YbZca hY XHJHJj' gY' YbU

**Version 2: Dec. 2017**